

Synd-EZ SS Kit Syndesmosis Repair

The Synd-EZ SS kit provides a solution for damage to the syndesmosis joint which connects the tibia and fibula bones of the lower leg. This simple yet effective device is designed to help stabilize the ankle and is intended for use with a stainless steel fracture plate.

Versatile stainless steel design

• Top hat is low profile and self centering on fracture plates

Adjustable knotless construct

• Simplifies surgical procedure

Counter traction suture & tension handle

 Integrated for predictably smooth tension of construct




Synd-EZ SS Kit Syndesmosis Repair

SYND-EZ SS KIT COMPONENTS

- 1. Top hat
- 2. White/blue adjustment suture
- 3. Solid blue counter-traction suture
- 4. Button
- 5. Solid white passing suture

- 6. White/black toggle suture
- 7. 3.5mm cannulated drill bit
- 8. 2.4mm passing pin
- 9. 1.5mm K-wire
- 10. Black tension handle





SYND-EZ SS ORDERING INFORMATION

Synd-EZ									
Part #	Width of But- ton (mm)	Length of Button (mm)	Description	Length of Loop (mm)	QTY/Box				
11224	3.2	12	Synd-EZ Stainless Steel Kit (SS)	-	1				
Synd-EZ Instrumentation									

Synd-EZ instrumentation									
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/ Reusable				
11272	-	Small Joint	Synd-Ez Drill Guide	Unsterile	Reusable				

Parcus Medical has joined Anika



Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL 1-941-755-7965 · www.parcusmedical.com ParcusCS@anika.com

