



Parcus Knee Instruments are designed to provide precise preparation of implant sites critical to performance and reproducible outcomes. Available individually or as part of the Parcus ACL Set.

## Features & Benefits

### Designed to meet the needs of various orthopedic procedures

- Arthroscopic and open techniques
- Suspensory and aperture fixation

### Laser depth marking (as applicable)

- Convenience in the operating room

ACL Set					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10077	-	Knee	Parcus, ACL Tray	Unsterile	Reusable

Tendon Harvesters					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10389	4.75	Knee	Closed tendor harvester	Unsterile	Reusable
10390	5.5	Knee	Closed tendor harvester	Unsterile	Reusable
10807	7	Knee	Closed tendor harvester	Unsterile	Reusable
10411	4.75	Knee	Slotted tendon harvester	Unsterile	Reusable
10421	5.5	Knee	Slotted tendor harvester	Unsterile	Reusable

Chondral Picks					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10554	-	Knee	Chondral pick (straight)	Unsterile	Reusable
10555	-	Knee	Chondral pick (15°)	Unsterile	Reusable
10556	-	Knee	Chondral pick (45°)	Unsterile	Reusable

Cannulated Headed Reamers					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10112	6	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10113	6.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10084	7	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10114	7.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10085	8	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10115	8.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10086	9	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10116	9.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10087	10	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10117	10.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10088	11	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10118	11.5	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable
10089	12	Knee	Headed reamer, w/2.5mm cannulation	Unsterile	Reusable



Cannulated Drill Bits					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10173	6	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10174	6.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10078	7	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10175	7.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10079	8	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10176	8.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10080	9	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10177	9.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10081	10	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10178	10.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10082	11	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10179	11.5	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable
10083	12	Knee	Drill bit, w/2.5mm cannulation	Unsterile	Reusable

Slotted Graft Sizers & Sizing Block					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10090	7	Knee	Slotted graft sizer	Unsterile	Reusable
10091	8	Knee	Slotted graft sizer	Unsterile	Reusable
10092	9	Knee	Slotted graft sizer	Unsterile	Reusable
10093	10	Knee	Slotted graft sizer	Unsterile	Reusable
10094	11	Knee	Slotted graft sizer	Unsterile	Reusable
10095	12	Knee	Slotted graft sizer	Unsterile	Reusable
10127	-	Knee	Graft sizing block		

Parcus Graft Station					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10107	-	Knee	Parcus graft station tray	Unsterile	Reusable
10110	-	Knee	Graft station fixed post	Unsterile	Reusable
10111	-	Knee	Graft station grasper and tensioner	Unsterile	Reusable
10143	-	Knee	Graft station base	Unsterile	Reusable
10162	-	Knee	Graft station fixed clamp	Unsterile	Reusable
20021	-	Knee	Graft station prep board	Unsterile	Reusable

Offset Guides, Tunnel Notcher, & Guide Pins					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10096	5	Knee	Femoral offset guide	Unsterile	Reusable
10097	6	Knee	Femoral offset guide	Unsterile	Reusable
10098	7	Knee	Femoral offset guide	Unsterile	Reusable
10838	4.3	Knee	Tunnel notcher	Unsterile	Reusable
10304	2.4	Knee	Drill tip guide pin, 250mm	Unsterile	Single-use
10301	2.4	Knee	Drill tip guide pin, w/eyelet, 250mm	Unsterile	Single-use
10307	2.4	Knee	Drill tip guide pin, 320mm	Unsterile	Single-use
10306	2.4	Knee	Drill tip guide pin, w/eyelet, 420mm	Unsterile	Single-use
10300	2.6	Knee	Trocar tip guide pin, w/eyelet, 380mm	Unsterile	Single-use



Tibial Plugs					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10512	13	Knee	Tibial plug, w/2.8mm cannulation	Unsterile	Reusable
10500	13	Knee	Tibial plug	Unsterile	Reusable

Ratcheting Handle Drivers					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10041	-	Knee	Hex titanium interference screw driver shaft, hudson, w/1.5mm cannulation	Unsterile	Reusable
10042	-	Knee	PEEK CF interference screw driver shaft, hudson, w/1.5mm cannulation	Unsterile	Reusable
20010	-	Knee	Ratcheting handle, w/3.8mm cannulation	Unsterile	Reusable

Fixed Handle Drivers					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10274	-	Knee	Hex titanium interference screw driver shaft, w/1.5mm cannulation	Unsterile	Reusable
10275	-	Knee	PEEK CF interference screw driver shaft, w/1.5mm cannulation	Unsterile	Reusable

Nitinol Guide Wires					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
20098	1.1	Knee	Nitinol guide wire, 400mm	Unsterile	Single-use
20015	1.5	Knee	Nitinol guide wire, 300mm	Unsterile	Single-use
20009	1.5	Knee	Nitinol guide wire, 400mm	Unsterile	Single-use

Parcus Pin Packs & Depth Gauge					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10057	-	Knee	GFS depth gauge	Unsterile	Reusable
10498	1.1	Knee	Parcus pin pack, includes: 10306, 10307, 20155 (1 per box)	Sterile	Single-use
10499	1.5	Knee	Parcus pin pack, includes: 10306, 10307, 20015 (1 per box)	Sterile	Single-use

Tibial Drill Guide					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10291	-	Knee	Cobra tooth drill guide aiming arm	Unsterile	Reusable
10299	-	Knee	Adjustable tibial drill guide	Unsterile	Reusable

