


PEEK CF Push-In Suture Anchors are recommended for knotted fixation in shoulder instability or foot and ankle.

Features & Benefits

Made of carbon fiber-reinforced PEEK-OPTIMA

- Non-absorbable, radiolucent, and MR safe
- Modulus of elasticity closely matching cortical bone



PEEK CF Push-In Suture Anchors

| Part # | Diameter (mm) | Length (mm) | Description | Material | QTY/Box |
|--------|---------------|-------------|--|----------|---------|
| 11339 | 3.5 | 8 | 3.5mm PEEK CF Push-In Suture Anchor w/ one 1.6mm Suture Tape (blu) | PEEK CF | 1 |
| 10323T | 3.5 | 8 | Push-in, suture anchor, w/ 2, #2 sutures (blu & wht/grn) | PEEK CF | 1 |
| 10407 | 3.5 | 8 | Push-in, suture anchor, w/ 2, #2 sutures (blu & wht/grn), w/needles (MO-6 short) | PEEK CF | 1 |

Instrumentation for PEEK CF Push-In Suture Anchors

| Part # | Diameter (mm) | Joint | Description | Sterile/Unsterile | Single-use/Reusable |
|--------|---------------|-------------|---|-------------------|---------------------|
| 10333U | 2.8 | Small Joint | Drill bit, w/positive stop | Unsterile | Reusable |
| 10538U | 2.8 | Shoulder | Drill bit w/positive stop | Unsterile | Reusable |
| 10482U | 3 | Small Joint | Drill bit w/positive stop | Unsterile | Reusable |
| 10537U | 3 | Shoulder | Drill bit w/positive stop | Unsterile | Reusable |
| 10823 | - | Small Joint | 8-Point slotted drill guide, w/beaver tail handle | Unsterile | Reusable |
| 10924 | 5.2 | Shoulder | 8-Point slotted knotless drill guide, w/beaver tail handle ¹ | Unsterile | Reusable |
| 10446 | - | Shoulder | 8-Point drill guide w/trocar tip obturator | Unsterile | Reusable |
| 10330 | - | Shoulder | V-Mouth drill guide w/trocar tip obturator | Unsterile | Reusable |
| 10406 | 3.4 | Shoulder | Blunt tip obturator ² | Unsterile | Reusable |
| 10771 | 3.3 | Shoulder | Conical tip obturator, w/1.8mm cannulation ² | Unsterile | Reusable |
| 10405 | 3.4 | Shoulder | Trocar tip obturator ² | Unsterile | Reusable |
| 10397 | 2.8 | Shoulder | Switching stick | Unsterile | Reusable |



Parcus Medical has joined Anika



Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 · www.parcusmedical.com

ParcusCS@anika.com

¹Accommodates 300mm drill guide. ²Use w/10446 & 10330